

FETAL MONITOR

JPD-300P

- 
Wired Connection
- 
Wireless Connection
- 
Alarm Function
- 
TFT Display Screen
- 
Clinical Event Marking Button



[10.1 inch]



[12.1 inch]



[7 inch]



3 in 1



Water Proof
IPX4

FEATURES

Standard Performance: FHR, TOCO, Fetal Movement, AUTO Fetal Movement.

Optional: Twins, Wireless transducer; Trolley; Wall Mount

Suitable: antepartum

- 90 degrees foldable TFT color screen
- 2 display modes: graphic and text modes
- 24 hours data storage
- 152mm thermal printing paper
- Compact and portable design

FETAL MONITOR


JPD-300P

ADVANTAGE

- Fetal Acoustic Stimulator (FAS) are built-in the FHR transducer, which is more convenient to stimulate sleeping fetus to speed up the test. (Patented).
- High sensitive and water-proof transducer
- CTG interpretation report
- Built-in battery lasting 4-hours
- Supports both cable and wireless network connection
- Connected to central monitoring system which up to 32 beds (Option)



Big font interface



Standard interface



Printer



Accessories



USB | Internet | Connector output



FHR1 | FHR2 | FM | TOCO



Power